

CELLS FOR PRECISION MEDICINE.

The epitope-independent Parsortix™ system enables isolation, identification and analysis of epithelial and mesenchymal CTCs

Anne-Sophie Pailhes-Jimenez, Senior Scientist

Molecular Diagnostics Europe
24 May 2018 - Lisbon



CELLS FOR PRECISION MEDICINE.




INTRODUCTION

Circulating Tumor Cells (CTCs)

Parsortix™ system

Introduction - Circulating Tumor Cells



Biology, detection, and clinical implications of circulating tumor cells. Joosse SA, Gorges TM, Pantel K. EMBO Mol Med. 2015 Jan.

- Reflective of metastatic tumor
- Source of biomarkers
- Non-invasive diagnostics
- Rare event: 1 CTC / 10⁸ blood cells


➡ **CTCs enrichment**

Parsortix™ technology




Platform technology for harvesting rare cells such as CTCs from blood


- based on **physical** cell properties (deformability)
- based on cell **size**
- **epitope-independent** technology



Parsortix cassette



Parsortix microfluidic system




Viabile cells captured by Parsortix system



CTCs are captured and enriched by the steps/gap in the Parsortix cassette

- Intact cell morphology



SKBR3 cells enriched by the Parsortix system and stained using Wright's and Giemsa staining

ANGLE's internal data

- Viable cells

→ cell culture
and *in vitro* testing

RESEARCH ARTICLE



Optimization and Evaluation of a Novel Size Based Circulating Tumor Cell Isolation System

Lei Xu^{1,2}, Xueying Mao¹, Ahmet Imrali¹, Ferriat Syed³, Katherine Mutsvangwa^{3,4}, Daniel Berney¹, Paul Cathcart^{4,5}, John Hines⁶, Jonathan Shamash^{3,4}, Yong-Jie Lu^{1*}

¹ Centre for Molecular Oncology, Barts Cancer Institute, Queen Mary University of London, London, United Kingdom, ² Department of Urology, Zhongshan Hospital, Fudan University, Shanghai, China, ³ Department of Medical Oncology, Barts Health NHS, London, United Kingdom, ⁴ Centre for Experimental Medicine, Barts Cancer Institute, Queen Mary University of London, London, United Kingdom, ⁵ Department of Urology, University College Hospital NHS, London, United Kingdom, ⁶ Department of Urology, Barts Health NHS, London, United Kingdom

* y.j.lu@qmul.ac.uk

Parsortix system – Recovery of cells




Pailhes-Jimenez et al. manuscript in preparation

➔ The Parsortix system is able to capture and harvest CTCs independently of cancer type and origin.

Parsortix system – Stability of cells




SKBR3 cell stability



The Parsortix system is compatible with different preservative blood collection tubes.

Parsortix system - Flexible analyses






CELLS FOR PRECISION MEDICINE.



DOWNSTREAM ANALYSIS

Some examples

Downstream analysis – Genomic analysis



- Molecular characterisation
- Specific to enriched subpopulations
- NGS/CGH/...
- Compare to ctDNA studies

In Situ Detection and Quantification of AR-V7, AR-FL, PSA, and KRAS Point Mutations in Circulating Tumor Cells
 Amin El-Heliebi,^{1,2*} Claudia Hille,^{2†} Navya Laxman,^{2†} Jessica Svedlund,^{2†} Christoph Haudum,^{1,4} Erkan Ercan,² Thomas Kroneis,² Shukun Chen,¹ Maria Smolle,⁶ Christopher Rossmann,² Tomasz Krzykowski,³ Annika Ahlford,^{3,6} Evangelia Darai,³ Gunhild von Amsberg,⁷ Winfried Alsdorf,⁷ Frank König,⁸ Matthias Löhr,⁹ Inge de Kruijff,¹⁰ Sabine Riethdorf,² Tobias M. Gorges,² Klaus Pantel,² Thomas Bauernhofer,^{4,5} Mats Nilsson,³ and Peter Sedlmayr¹



A Novel Workflow to Enrich and Isolate Patient-Matched EpCAM^{high} and EpCAM^{low/negative} CTCs Enables the Comparative Characterization of the PIK3CA Status in Metastatic Breast Cancer
 Rita Lampignano, Liven Yang, Martin H. D. Neumann, André Franken, Tanja Fehn, Dieter Niederacher and Hans Neubauer*



Clinical Chemistry 62:11
000-000 (2016)

Cancer Diagnostics


Accession of Tumor Heterogeneity by Multiplex Transcriptome Profiling of Single Circulating Tumor Cells

Tobias M. Gorges,^{1†} Andra Kuske,^{1†} Katharina Röck,¹ Oliver Mauermann,¹ Volkmar Müller,² Sven Peine,² Karl Verpoort,⁴ Vendula Novosadova,⁵ Mikael Kubista,^{5,6} Sabine Riethdorf,¹ and Klaus Pantel^{1†}

A novel microfluidic platform for size and deformability based separation and the subsequent molecular characterization of viable circulating tumor cells

G.E. Hvhichia¹, Z. Parveen¹, C. Wagner¹, M. Janning^{2,3,4}, J. Quidde^{3,4}, A. Stein^{3,4}, V. Müller⁵, S. Loges^{2,3,4}, R.P.L. Neves⁶, N.H. Stoecklein⁶, H. Wikman², S. Riethdorf², K. Pantel² and T.M. Gorges²

Downstream analysis - FISH



➤ Genomic alterations

Personalized Medicine and Imaging | **Clinical Cancer Research**


The Novel Association of Circulating Tumor Cells and Circulating Megakaryocytes with Prostate Cancer Prognosis

Lei Xu^{1,2}, Xueying Mao¹, Tianyu Guo¹, Pui Ying Chan³, Greg Shaw⁴, John Hines⁴, Elzbieta Stankiewicz¹, Yujin Wang¹, R. Tim D. Oliver¹, Amar Sabri Ahmad², Daniel Berney¹, Jonathan Shamash³, and Yong-Jie Lu¹

Clinical Chemistry 64:3 000-000 (2018)


In Situ Detection and Quantification of AR-V7, AR-FL, PSA, and KRAS Point Mutations in Circulating Tumor Cells

Amin El-Heliebi,^{1†} Claudia Hille,^{2†} Navya Laxman,^{3†} Jessica Svedlund,^{3†} Christoph Haudum,^{1,4} Erkan Ercan,¹ Thomas Kronels,¹ Shukun Chen,³ Maria Smolle,⁴ Christopher Rossmann,⁵ Tomasz Krzykowski,⁴ Annika Ahlford,^{3,6} Evangelia Darai,³ Gunhild von Amsberg,⁷ Winfried Aisdorf,⁷ Frank König,⁸ Matthias Lohr,⁹ Inge de Kruijff,¹⁰ Sabine Riethdorf,² Tobias M. Gorges,² Klaus Pantel,² Thomas Bauernhofer,¹¹ Mats Nilsson,¹² and Peter Sedlmayr¹




SKBR3 cells enriched by the Parsortix system and hybridized using PathVysion Her-2 DNA Probe Kit
 ANGLE's internal data

Downstream analysis - Immunofluorescence



- Morphology
- Real-time visualisation of expressed protein
- Molecular characterisation
- **DAPI⁺CK⁺CD45⁻**



SKBR3 cells enriched by the Parsortix system and stained using CK/EpCAM Ab.
Pailhes-Jimenez et al. manuscript in preparation

RESEARCH ARTICLE
Optimization and Evaluation of a Novel Size Based Circulating Tumor Cell Isolation System

Lei Xu^{1,2}, Xueying Mao¹, Ahmet Imrali¹, Ferriat Syed³, Katherine Mutsavangwa^{3,4}, Daniel Berney¹, Paul Cathcart^{4,5}, John Hines⁶, Jonathan Shamash^{2,4}, Yong-Jie Lu^{1*}

Analyst, 2016 Jan 21;141(2):669-78. doi: 10.1039/c5an02156a.

Clinical evaluation of a novel microfluidic device for epitope-independent enrichment of circulating tumour cells in patients with small cell lung cancer.

Chudziak J¹, Burt DJ¹, Mohan S¹, Rothwell DG¹, Mesquita B¹, Antonello J¹, Dalby S¹, Ayub M¹, Priest L¹, Carter L¹, Krebs MG¹, Blackhall F², Dive C¹, Brady G¹.

www.impactjournals.com/oncotarget/ Oncotarget, Advance Publications 2017


Comparison of isolation platforms for detection of circulating renal cell carcinoma cells

Yvonne Maertens¹, Verena Humberg¹, Franziska Erlmeier², Sandra Steffens¹, Julie Steinestel¹, Martin Bögemann¹, Andres Jan Schrader¹ and Christof Bernemann¹



Article
A Novel Workflow to Enrich and Isolate Patient-Matched EpCAM^{high} and EpCAM^{low/negative} CTCs Enables the Comparative Characterization of the PIK3CA Status in Metastatic Breast Cancer


Rita Lampignano, Liwen Yang, Martin H. D. Neumann, André Franken, Tanja Fehm, Dieter Niederacher and Hans Neubauer*




Development of a Parsortix workflow: In-cassette staining



- Automated method
- User-friendly
- Uses a slide-sized cassette
- Minimal cell loss (*spiking experiment*)



Development of a Parsortix workflow: In-cassette staining



*Cancer cells captured by the Parsortix system and stained using CK/EpCAM Abs.
Pailhes-Jimenez et al. manuscript in preparation*


Parsortix in-cassette staining Patient samples




12 healthy donor blood samples
30 **stage I-III BC** patient samples

Separation and staining using **Parsortix** instrument

CTC detection:
0% healthy donor
63% **stage I-III BC** samples



➔ CTCs are detectable in **63% stage I-III BC** patient samples.




Cancer cells captured by the Parsortix system and stained using CK/EpCAM Abs.
Pailhes-Jimenez et al. manuscript in preparation

CTC subpopulations detected by Parsortix in-cassette staining



Example: Detection of EMT-associated markers



Positive selections:


- Epithelial CTCs (CK/EpCAM based)= green
- Mesenchymal CTCs (CONFIDENTIAL)= orange

Negative selection:

- White Blood Cells (CD45 detection)= red

SKBR3 and Hs578T cells spiked into healthy blood, separated and stained using Parsortix technology.
Pailhes-Jimenez et al. manuscript in preparation

Example: Detection of phosphorylated proteins



Positive selections:

- Epithelial CTCs (CK/EpCAM based)= green
- Phosphorylated ERK = orange

Negative selection:

- White Blood Cells (CD45 detection)= red

PDBu-treated A549 cells spiked into healthy blood, separated and stained using Parsortix technology.
Pailhes-Jimenez et al. manuscript in preparation

Parsortix in-cassette staining workflow=
semi-automated and flexible method allowing for the isolation and
identification of CTCs by immunofluorescence

➤ Semi-automated method:

- In-cassette staining using Parsortix instrument
- Fast identification of captured cells (few hours)
- Uses a standard fluorescent microscope

➤ Flexible method:

- Use of different blood stability tubes
- Combination of targets/biomarkers (6 reagent tubes available)
- Possibility to harvest cells after in-cassette staining

Stained cells can be harvested and used for further analysis

Accession of Tumor Heterogeneity by Multiplex Transcriptome Profiling of Single Circulating Tumor Cells

Tobias M. Gorges,^{1,†} Andra Kuske,^{1,†} Katharina Röck,¹ Oliver Mauermann,¹ Volkmar Müller,² Sven Peine,³ Karl Verpoort,⁴ Vendula Novosadova,⁵ Mikael Kubista,^{5,6} Sabine Riethdorf,¹ and Klaus Pantel^{1*}

➤ Single cell multiplexing RNA profiling

In Situ Detection and Quantification of AR-V7, AR-FL, PSA, and KRAS Point Mutations in Circulating Tumor Cells

Amin El-Heliebi,^{1,†} Claudia Hille,^{2†} Navya Laxman,^{3†} Jessica Svedlund,^{3†} Christoph Haudum,^{1,4} Erkan Ercan,¹ Thomas Kroneis,¹ Shukun Chen,¹ Maria Smolle,⁴ Christopher Rossmann,⁵ Tomasz Krzywkowski,³ Annika Ahlford,^{3,6} Evangelia Darai,³ Gunhild von Amsberg,⁷ Winfried Alsdorf,⁷ Frank König,⁸ Matthias Lühr,⁹ Inge de Kruijff,¹⁰ Sabine Riethdorf,² Tobias M. Gorges,² Klaus Pantel,² Thomas Bauernhofer,^{4,5} Mats Nilsson,³ and Peter Sedlmayr¹

➤ *In-situ* mutation detections (AR-V7)

Personalized Medicine and Imaging Clinical Cancer Research

The Novel Association of Circulating Tumor Cells and Circulating Megakaryocytes with Prostate Cancer Prognosis

Lei Xu^{1,2}, Xueying Mao¹, Tianyu Guo¹, Pui Ying Chan³, Greg Shaw⁴, John Hines⁴, Elzbieta Stankiewicz², Yuqin Wang¹, R. Tim D. Oliver¹, Amar Sabri Ahmad⁵, Daniel Berney¹, Jonathan Shamash³, and Yong-Jie Lu¹

➤ FISH

International Journal of Molecular Sciences MDPI

Article

A Novel Workflow to Enrich and Isolate Patient-Matched EpCAM^{high} and EpCAM^{low/negative} CTCs Enables the Comparative Characterization of the PIK3CA Status in Metastatic Breast Cancer

Rita Lampignano, Liwen Yang, Martin H. D. Neumann, André Franken, Tanja Fehm, Dieter Niederacher and Hans Neubauer *

➤ Single cell + sequencing




CELLS FOR PRECISION MEDICINE.



CONCLUSIONS

Conclusions



Thank you



ANGLE's R&D team and selected collaborators



Parsortix™ system enables isolation, identification and analysis of CTCs



Questions?



ANGLE Europe Ltd
10 Nugent Road
The Surrey Research Park
Guildford GU2 7AF
United Kingdom

www.ANGLEplc.com

ANGLE North America Inc
3711 Market Street
University Science Center 8th floor
Philadelphia PA 19104
USA

ANGLE Biosciences Inc
50 Ronson Drive, Suite 105
Toronto
Ontario M9W 1B3
Canada