


Exploring the dual analysis of circulating tumour cells and circulating tumour DNA



Anne-Sophie Pailhes-Jimenez
Senior Director, R&D
ANGLE Europe Ltd



Disclaimer




I am an employee from ANGLE Europe Limited.

The opinions expressed during this presentation are my own and ANGLE's and may not represent the opinions of Illumina.

Any uses of ANGLE and Illumina products described in this presentation may be uses that have not been cleared or approved by the FDA. or any other applicable regulatory body.

Illumina is compensating me to speak at this event.


The potential of liquid biopsy throughout the patient care pathway



Liquid Biopsy and CTCs


Blood as liquid Biopsy:

- **Minimally-invasive**
- **Repeatable**
- **Representative** of cancer heterogeneity
- **Versatile**
- **Informative**



Klaus Pantel and Catherine Alix-Panabieres.
 Liquid biopsy and minimal residual disease — latest advances and implications for cure.
 Nat Rev Clin Oncol **16**, 409–424 (2019). <https://doi.org/10.1038/s41571-019-0187-3>

CTCs and ctDNA provide complementary information



Hypothesis:

Dual analysis provides critical information not otherwise available on cancer progression and for targeted treatment selection



“**CTCs**, as living cells that are active in the metastatic process, can provide **prospective insight into a patient’s cancer**.”

In comparison, **ctDNA** derived mainly from dead and dying cells, provides important but **historical** information on patient disease.”

Prof. Evi Lianidou

Head of the Molecular Diagnostics Laboratory
National and Kapodistrian University of Athens

Circulating tumour cells (CTCs)



Biomarkers circulating in the bloodstream

CTC and CTC clusters:

- Reflective of **metastatic** tumor
- Source of **biomarkers** (suitable for multiomics analysis)
- **Rare** event: 1 CTC / 10^9 blood cells

➔ Importance of **upstream workflow** to generate quality sample.

Parsortix[®] PCI system



The first FDA cleared / CE-IVD medical device for the capture and harvest of circulating tumor cells from metastatic breast cancer patient blood for user-validated subsequent analysis.



Intended use

“The Parsortix[®] PCI system is an in vitro diagnostic device intended to enrich circulating tumor cells (CTCs) from peripheral blood collected in K₂EDTA tubes from patients diagnosed with metastatic breast cancer. The system employs a microfluidic chamber (a Parsortix cell separation cassette) to capture cells of a certain size and deformability from the population of cells present in blood. The cells retained in the cassette are harvested by the Parsortix PCI system for use in subsequent downstream assays.”


The end user is responsible for the validation of any downstream assay. The standalone device, as indicated, does not identify, enumerate or characterize CTCs and cannot be used to make any diagnostic/prognostic claims for CTCs, including monitoring indications or as an aid in any disease management and/or treatment decisions.




*FDA Clearance for use with metastatic breast cancer patients only

Parsortix systems: capturing and harvesting living cancer cells

Plan view



Patented multifold and separation step




Platform technology

The Parsortix system harvests **cancer cells from blood** based on their larger size and lack of deformability

Other cells can be captured:

- **white blood cells** associated with the tumor microenvironment
- **megakaryocytes** (frequency may relate to cancer)
- **fetal cells** from maternal blood



CTC and WBC cluster

Key


- Captured CTCs
- White blood cells
- Red blood cells
- Blood flow

Workflow



Workflow allowing for CTCs and ctDNA analysis

Dual analysis of biomarkers present in plasma and cellular component



<https://usc norriscancer.usc.edu/liquid-biopsy-research-core/>

Proof Of Concept work:

No CTC loss following plasma removal, prior to Parsortix enrichment.




Figure 5. CTC detection in MBC patients blood samples processed with or without plasma. Dot plot shows mean ± SEM of the number of A) CTCs and B) CTC clusters identified in each donor per condition. No statistically significant difference was observed, $ns = p \geq 0.05$, Paired T-test. C) Table showing number of donors included in each cohort (N), percentage (%) of donors with ≥1 CTC, range, mean and median of CTCs captured within the positive group, percentage of donors with ≥1 CTC cluster, range of CTC clusters per donor, and range of number of CTCs per cluster; D) Representative images of a single mesenchymal CTC (top) and a cluster of mesenchymal CTCs (bottom). Epithelial markers (FITC) in green, Mesenchymal markers (Cy7) in magenta, Blood lineage markers (Cy5) in red or white in the merge, Nucleus (DAPI) in blue.


<https://angleplc.com/wp-content/uploads/2023/05/PTX-Poster-A-ISLB-2022-Bidetect-002.pdf>



Can we allow for dual detection of targets/variants which may have an impact on patient care?

Molecular analysis of CTCs and ctDNA

ctDNA and Parsortix-derived CTC-DNA can be analysed from single blood sample across multiple cancer types to detect mutations using **Illumina Cell-Free DNA prep with enrichment combined with custom panel targeting Lung Cancer.**



Illumina Cell-Free DNA prep with Enrichment and custom 79 gene panel enables a flexibility in the targeted genes.

DNA dual analysis

Assay performance (1/2)

Study design for evaluation of performance of DNA dual analysis

Contrived samples were used to assess the analytical performance of the assay:


- Two cell lines with known mutations, were spiked into Parsortix-enriched samples
- DNA samples from CTC were extracted and processed following Illumina Cell-Free DNA prep with Enrichment procedure.

Table: expected mutations to be detected in spiked samples

Cell line 1	Cell line 2
pTEN p.T167A	EGFR p.T790M
BRAF p.V600E	EGFR p.L858R
CDK4 p.R24C	PIK3CA p.G118D
	TP53 p.R273H

DNA dual analysis

Assay performance (1/2)



DNA dual analysis of Patients with Lung cancer

Patient data 1/4



2 cohorts of Lung Cancer patient

- **1st cohort: 8 Lung Cancer patient samples (no treatment)**
- **2nd cohort: 19 Lung Cancer patient samples (may be under treatment)**

✓ **Good quality** samples

✓ Overall, **93% of samples with at least one oncogenic gene detected.**

DNA dual analysis of Patients with Lung cancer

Patient data 2/4

1st cohort - Lung cancer panel (n=8):

CTC(average)= 49 and CTC(median)= 20



100% patients showed oncogenic mutations signal for one or more targets.

32 genes (74 variants) seem to be specifically detected in this cohort:

Table: example of detected mutations in specific analyte

CTC	cfDNA
CHEK1*	MSH6
FBXW7*	POLE
MLH1	BRAF
ESR1	APC
ROSI	ATM

* predominant



Dual analysis of Patients with Lung cancer

Patient data 3/4

2nd cohort - Lung cancer panel (n=19):

CTC(average)= 4 and CTC(median)= 3

90% of patients with at least one oncogenic mutations in CTC for one or more targets.

38 genes (214 variants) have been specifically detected in this cohort:





Table: example of detected mutations in specific analyte

CTC	cfDNA
EpCAM	CDKN2A
MYCN	NTRK
FOXL2	PIK3CA
HRAS	APC
CHEK1*	RBI*

* predominant


DNA dual analysis of Patients with Lung cancer

Patient data 4/4



Full cohort - Lung cancer panel (n=27):

43 genes (247 variants) have been specifically detected in this cohort:



Analysing both analytes give access to **twice** as many **mutations**, which could be investigated and/or use as “potential” biomarkers.

DNA dual analysis of Patients with Lung cancer

Validation of data using another technology

➔ **Validation** of patient data using another cohort of patient samples (under treatment – n=10)

Use of a pan-CANCER NGS panel, with blocker displacement amplification

20 genes (27 variants) have been specifically detected in this cohort:





Table: example of detected mutations in CTC fractions

CTC	cfDNA
PiK3CA	MYC
MTOR	MET
ESR I	EGFR*

** predominant*




Confirmation of the utility of using CTC as well as cfDNA to have access to a larger number of potential biomarkers

Conclusions



Parsortix DNA dual analysis strategy

CONCLUSION



- ✓ Successful detection of potential actionable targets in cfDNA and CTC, from the same sample.
- ✓ CTC samples provide a **unique detection** of actionable targets in most of the patient samples.
- ✓ This study highlights the potential for multiomic analysis.
- ✓ More data to come...

Parsortix system is allowing for CTC enrichment and **DNA analysis from plasma and CTC, from the same blood sample, with high analytical sensitivity and specificity.**

Dual analysis provides the **possibility to access (new) biomarkers**, which may help guide patient care.

Thank you for your attention!

ANGLE plc

10 Nugent Road

Surrey Research Park

Guildford GU2 7AF

United Kingdom

www.angleplc.com

Anne-Sophie PAILHES-JIMENEZ

a.pailhes-jimenez@angleplc.com

